

we deliver
you discover



THE CAPMATE PLUS

single tube capping device

*A consistently
secure device that
automatically seals
your tubes and vials
with screw caps.*

ELIMINATES
FATIGUE IN
THE LAB!

Features

- Caps 20-30 tubes per minute
- Bowl feeder automatically moves caps into capping head
- Digital cap counter tracks how many tubes have been capped
- Automatic button sends tubes into unit, capping them automatically

benefits

- Expedites and automates a segment of laboratory workflow
- Eliminates tedious and non-essential tasks from the lab
- Increases productivity and innovation
- Provides more time to focus on complicated scientific problems
- Reduces operator hand fatigue and prevents occupational hazards like carpal tunnel syndrome
- Small footprint

specifications

- Unit size: 14.1 x 8.7 x 19.7 in
- Weight: 19.4 lbs
- Power supply: 115 V

