

we deliver
you discover



INTRODUCING THE TBS-1080!

The only chlorine dioxide gas generation system that meets industry standards AND fits inside a BSC.



FEATURES AND BENEFITS

- Small, lightweight & portable (only 12 lbs)
- User-friendly operation
- Large LCD screen shows progress of disinfection
- Disinfects in 90 minutes
- 2 steps prevent accidental start
- Fits inside a biological safety cabinet
- Effective on all forms of microorganisms, viruses and spores
- Uses prepackaged reagents and refillable chemical cartridges
- Chlorine dioxide gas recognized by EPA as acceptable disinfection method
- No chlorine gas generated

FITS
INSIDE A
BIOSAFETY
CABINET!