

SCREW CAP RECAPPER CS500



The Micronic Screw Cap Recapper CS500 is developed in conjunction with several leading laboratories. The handheld device offers users the ability to cap, decap or recap a row of 8 tubes in less than 5 seconds, improving productivity!

The Micronic Screw Cap Recapper CS500 requires no battery charging and is always available for use. Whether you are looking to recap tubes with the same or new caps it is easy to position the screw caps above tubes using the CS500. When removing screw caps, the large ejection handle makes discarding the old caps a simple single action task.

Tube compatibility: The Micronic Recapper CS500 is compatible with Micronic standard Screw Caps (internal and external thread) or Low Profile Screw Caps.

SUPERIOR FEATURES

- Highly durable solution to cap, decap or recap a row of 8 tubes at once
- Each screw cap is reproducibly tightened securely and uniformly
- Caps a row of 8 tubes in 5 seconds and a whole rack of 96 tubes in only 1 minute - 10 times faster than by hand!
- Minimizes contamination
- Easy positioning of screw caps above the tubes
- Ergonomic lightweight solution for left and right handed use
- Caps can be easily discarded using the large ejection handle
- Very small footprint for bench top use
- No need for battery recharging
- Models available for Micronic 96-format, 48-format and 24-format Screw Caps, as well as caps from other brands
- 3-year warranty




Highly durable solution



(DE) CAPPERS

SCREW CAPS

SPECIFICATIONS



Compatibility:	All ANSI / SLAS standard 96-format racks
Cycle time:	Less than 5 seconds per row
Dimensions:	17.8 x 16.1 x 71.0 cm / 7.01 x 6.34 x 27.95 inches
Weight:	5.5 kg / 12.0 lbs
Voltage:	115 - 230 V

ORDERING INFORMATION

Catalog number	Description	Unit
Screw Cap Recappers		
MP35420	Micronic Screw Cap Recapper CS500 for Micronic Screw Caps (TRE)	1 unit
MP35421	Micronic Screw Cap Recapper CS500 for Micronic Screw Caps (TRI)	1 unit
MP35422	Micronic Screw Cap Recapper CS500 for Micronic Low Profile Screw Caps (TRI)	1 unit
MP35423	Micronic Screw Cap Recapper CS500 for Matrix Screw Caps (TRI)	1 unit
MP35424	Micronic Screw Cap Recapper CS500 for FluidX Screw Caps (TRE)	1 unit
MP35425	Micronic Screw Cap Recapper CS500 for FluidX Screw Caps (TRI)	1 unit
MP35426	Micronic Screw Cap Recapper CS500 for FluidX 2K Screw Caps (TRI)	1 unit
MP35427	Micronic Screw Cap Recapper CS500 for LVL Screw Caps (TRE)	1 unit
MP35428	Micronic Screw Cap Recapper CS500 for LVL Screw Caps (TRI)	1 unit
MP35441	Micronic Screw Cap Recapper CS510 for Micronic Screw Caps (TRI)	1 unit
MP35446	Micronic Screw Cap Recapper CS520 for Micronic Screw Caps (TRI)	1 unit
MP51500	Micronic disposable cap carrier for 96-well format tubes (empty)	10 carriers
MP51501	Micronic disposable cap carrier for 48-well format tubes (empty)	10 carriers
MP51502	Micronic disposable cap carrier for 24-well format tubes (empty)	10 carriers

AWARD WINNING DESIGNS

The Micronic recappers CS500 and CS700 have been awarded with a good industrial design award (GIO Award)! An independent professional jury gave this award to both the CS500 and the CS700 during Dutch Design Week 2017. This means that the designs are original, functional, ergonomic and respectful towards the environment. Moreover, the GIO award indicates that the new designs add value to end-users. The CS500 and CS700 are the first two products in the line of the Micronic equipment with a new design. More newly designed equipment is coming soon, so keep in touch!



2017:
good industrial
design