

MiniVortex Mixer

◆ Description

The mini vortex mixer is the ideal vortex mixer for every lab with its compact size, ease of use, and its economic price point. The 3,000 rpm fixed speed ensures there is no sacrificing any mixing power.

◆ Features

- Extremely compact and easy operation.
- Head is made of durable material to ensure long life.
- Refined steel base is designed to ensure product stability.
- Powerful vortexing for tubes up to 30mm.
- Choice of different colored working head.



Mini Vortex Mixer
VM-03UA



Mini Vortex Mixer
VM-03RUA



Mini Vortex Mixer
VM-03GUA






Mini Vortex Mixer
VM-03PUA

◆ Specifications

Name	Mini Vortex Mixer			
Model	VM-03UA	VM-03RUA	VM-03GUA	VM-03PUA
Color	Blue	Red	Green	Purple
For regions other than North America, check the model number in the "Order Information Table"				
Speed	3000rpm			
Speed Accuracy	± 10%			
Oscillation Mode	Circumference Oscillation			
Orbital Diameter	Ø4.5mm (Ø0.18")			
Capacity	Max. 80ml			
Accessories included in Standard Packaging	NA			
Tube Diameter	≤Ø30mm (Ø1.18")			
Ambient Temperature	5~40°C			
Relative Humidity	≤80%			
Power Supply	AC100~240V, 50 / 60Hz			
Max. Power	5W			
Dimensions (ΦxH)	96x85mm (3.78x3.35")			
Net. Weight	0.55kg (1.21 lb)			
Protection Level	IP21			
Certifications	CE, RoHS, WEEE			

◆ Order Information Table

Order Model	Power Supply	Plug Type
VM-03UA	AC100-240V, 50 / 60Hz	US 
VM-03RUA		
VM-03GUA		
VM-03PUA		
VM-03UW		Eur, UK, AU, US 
VM-03RUW		
VM-03GUW		
VM-03PUW		
VM-03USA		
VM-03RUSA		
VM-03GUSA	SA 	
VM-03PUSA		