



Pipette Tips

For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.

Universal Pipette Tips with Advanced Features and Premium Materials



PT-10-R

10 μ L Universal Tip


Material: PP (Polypropylene)
Length: 32mm
Capacity: 10 μ L
Order QTY: 10 Racks per pack
(96 tips in each rack)



PT-200-R

200 μ L Universal Tip

Material: PP (Polypropylene)
Length: 51mm
Capacity: 200 μ L
Order QTY: 10 Racks per pack
(96 tips in each rack)



PT-1000-R

1000 μ L Universal Tip

Material: PP (Polypropylene)
Length: 78mm
Capacity: 1000 μ L
Order QTY: 10 Racks per pack
(100 tips in each rack)

Product Features:

- Produced in a Class 100,000 clean room
- Pharmaceutical grade Polypropylene
- DNase-free, RNase-free, pyrogen-free
- No endotoxin or PCR inhibitor
- Works with various pipettes such as Gilson, Eppendorf, Thermo (Finn)
- Hydrophobic filter avoids cross contamination
- Low retention tips reduce residue on the inner wall to ensure accurate results
- Gamma Irradiated Sterile

