

we deliver
you discover



THREE
LEVELS OF
SERVICE
AVAILABLE!

SERVICE PLANS FOR CS700-SERIES SCREW CAP RECAPPER

Four plans available to fit your needs.

1A



STANDARD PM SERVICE PLAN

Ship unit to workshop

- One (1) scheduled PM inspection (recommended every 12 months)
- Performance test
- Thorough cleaning and lubrication of machine axis, spindles and guidings
- Functionality testing of machine axis, motor, push-off unit, optional tube detection sensor, slip couplings, bits and o-rings
- Inspection of screw head, fuses, guides, o-rings, springs, power cords and plugs (adjustment, replacement or repair of these parts done as necessary)
- Exclusions: all repair parts

1B




STANDARD PM SERVICE PLAN

Repair unit on-site

- One (1) scheduled PM inspection (recommended every 12 months)
- Travel expenses related to PM visit
- Performance test
- Thorough cleaning and lubrication of machine axis, spindles and guidings
- Functionality testing of machine axis, motor, push-off unit, optional tube detection sensor, slip couplings, bits and o-rings
- Inspection of screw head, fuses, guides, o-rings, springs, power cords and plugs (adjustment, replacement or repair of these parts done as necessary)
- Exclusions: all repair parts

2



LEVEL 2 SERVICE CONTRACT

Repair unit on-site

- One (1) annual PM inspection (includes all offerings in Standard PM)
- One emergency, on-site service visit
- Service visit includes repair parts*, labor cost, and travel
- *Exclusions: screw head

3



LEVEL 3 SERVICE CONTRACT

Repair unit on-site

- Includes all offerings in Level Two
- Extra emergency on-site service visit
- Micronic software upgrades
- Service visit includes repair parts*, labor cost and travel
- *Exclusions: screw head

Need a quote? Have a question? Contact us today:



724-941-6411



info@nbsscientific.com



www.nbsscientific.com