

ExpellPlus Filter Tips

COMPATIBILITY CHART

For tips compatibility with other pipette brands, please contact info@capp.dk



5030030C
5131030C



5030060C
5131060C



5030062C
5131062C



5030066C
5131066C



5030090C
5131090C



5130150C
5131150C
5130120C

Capp Bravo, ecopipette and
CappAero mechanical pipette

0.2-2µl
0.5-10µl

0.2-2µl
0.5-10µl

2-20µl

5-50µl
10-100µl

2-20µl
5-50µl
10-100µl
20-200µl

100-1000µl

Capp electronic pipettes

x

2-20µl

x

x

10-200µl

100-1000µl

Biohit eLine, mLine, Proline Plus

0.1-3µl
0.5-10µl

0.1-3µl
0.5-10µl

2-20µl

2-20µl
10-100µl

10-100µl
20-200µl

100-1000µl

Brand Transferpette, Transferpette S,
Transferpette Electronic

0.5-10 µl

0.5-10µl

2-20µl

x

20-200µl

100-1000µl

Eppendorf Reference, Research Plus

0.5-10µl

0.1-2.5µl
0.5-10µl

2-20µl

x

20-200µl

100-1000µl

Finnpipette Digital

0.5-10 µl

0.5-10µl

2-20 µl

5-50µl

5-50µl
20-200µl

100-1000µl

Gilson Pipetman

0.2-2µl
1-10µl

0.2-2µl
1-10µl

2-20µl

2-20µl
10-100µl

2-20µl
10-100µl
20-200µl

100-1000µl

Rainin Classic, universal shaft (non LTS)

0.1-2µl
0.5-10µl

0.1-2µl
0.5-10µl

2-20 µl

2-20µl
10-100µl

2-20µl
10-100µl
20-200µl

100-1000µl

VWR, Labnet, HTL mechanical pipettes

0.1-2µl
0.5-10µl

0.1-2µl
0.5-10µl

2-20 µl

2-20µl
10-100µl

2-20µl
10-100µl
20-200µl

100-1000µl

Expell & ExpellPlus Tips

COMPATIBILITY CHART

For tips compatibility with other pipette brands, please contact info@capp.dk

					
	5130010C 5030010C 5030015C 5030020C 5030035C	5030040C 5030045C 5030050C	5130070C 4130075C 5030070C 5030075C 5030080C 5030095C	5130100C 5130105C 5130110C	5130130C 4130135C 5130135C 5130140C 5130155C
Capp Bravo, ecopipette and CappAero mechanical pipette	0.2-2µl 0.5-10µl	0.2-2µl 0.5-10µl	2-20µl 5-50µl 10-100µl 20 - 200µl	2-20µl 5-50µl 10-100µl 20-200µl 30-300µl	100-1000µl
Capp electronic pipettes	2-20µl	2-20µl	10-200µl	10-200µl	100-1000µl 100-1200µl
Biohit eLine, mLine, Proline Plus	0.1-3µl 0.5-10µl	0.1-3µl 0.5-10µl	2-20µl 10-100µl 20 - 200µl	2-20µl 10-100µl 20 - 200µl	100-1000µl
Brand Transferpette, Transferpette S, Transferpette Electronic	0.5-10µl	0.5-10µl	2-20µl 10-100µl 20 - 200µl	20-200µl 15-300µl	100-1000µl
Eppendorf Reference, Research Plus	0.1-2.5µl 0.5-10µl	0.1-2.5µl 0.5-10µl	2-20µl 10-100µl 20 - 200µl	10-100µl 20-200µl 30-300µl	100-1000µl
Finnpipette Digital	0.5-10µl	0.5-10µl	2-20µl 5-50µl 20-200µl	5-50µl 20-200µl 30-300µl	100-1000µl 200-1000µl
Gilson Pipetman	0.2-2µl 1-10µl	0.2-2µl 1-10µl	2-20µl 20-100µl 50-200µl	2-20µl 20-100µl 50-200µl	200-1000µl
Rainin Classic, universal shaft (non LTS)	0.1-2µl 0.5-10µl	0.2-2µl 1-10µl	2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	100-1000µl
VWR, Labnet, HTL mechanical pipettes	0.1-2µl 0.5-10µl	0.2-2µl 1-10µl	2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	100-1000µl

INFO@CAPP.DK

WWW.CAPP.DK